MDBio Inc.

產品手冊

Product information sheet

免费服务专线:4000-855-868 网站:www.mdbio.com.cn 电子邮件:mdbio@mdbio.com.cn

Product Name & Product No.

Product Name: X-GlcA 【129541-41-9】 Product No.: X003

Spec.

Enzymatic Colorimetric Reagent, For GUS Detection 5-Bromo-4-chloro-3-indolyl-β-D-Glucuronic Acid, Sodium Salt C14H12BrClNO7Na M.W.: 444.6 Assay: 98.0% min. Solubility(2% Soln. in DMF): Clear and Colorless

Description

X-GlcA, 5-Bromo-4-chloro-3-indolyl-b-D-Glucuronic acid is a substrate for β -D-Glucuronidase (GUS) encoded by the gus A gene. The substrate is used as a qualitative histochemical marker of specific GUS expression in cells and tissue. X-GlcA is cleaved by GUS at the β 1 glucuronic bond between glucuronic acid and the 5-Bromo-4chloro-3-indolyl part of X-GlcA via hydrolysis. The cleavage of X-GlcA results in the precipitation of water insoluble blue dichlorodibromo-indigo precipitate at the site of enzymatic cleavage. Colour formation requires three separate reactions. After enzymatic turnover, the released indoxyl derivative dimerizes and is subsequently oxidized to the final indigo dye.

Preparation Instruction

1% (w/v) X-GlcA, sodium salt in sterile N, Ndimethylformamide or deionized water (10 mg/mL).

Note:

X-GlcA, sodium salt is unstable in water. Use immediately after dissolution. X-GlcA solution should be added to media after it is cooled below 55 ∞ .

Applications

 Suitability for uidA (GUS) Detection in *Escherichia coli* uidA⁻ and uidA⁺ *Escherichia coli* cells were
streaked on separate Luria plates containing 1%
(w/v) peptone , 0.5% (w/v) yeast extract, 0.5%
(w/v) NaCl , 0.7% (w/v) agar and 0.33 mg/mL
X-GlcA, sodium salt. The plates were incubated
overnight at 37 °C. After 24-48 hours at 37 °C,
the uidA⁺ cells produced blue colonies
indicating the expression of the β-glucuronidase
gene and the uidA⁻ cells produced white
colonies indicating the absence of expression.

2. Suitability as Substrate for β -Glucuronidase A qualitative enzymatic assay was run in 0.1 M acetate buffer, pH 5.0 at 37 °C and 0.1 mg/mL X-GlcA with 6,000 units of β -glucuronidase in a 3 mL reaction mix. Suitability is determined by measuring the change in absorbance at 400 nm in one minute. A positive rate represents suitability.

Storage

Store at 4° C; Protect from light and humidity. For long term storage hold at -20°C.